

KEM

EN

TEC

Nordic

Ultrafast – Simply Brilliant




InstantBlue™ Ultrafast protein stain

Remove gel




15 minutes

Add stain



Finished



Ultra Quick

Results within
15 minutes

Easy

One step protocol

- No washing
- No fixing
- No heating
- No de-staining

Sensitive

≥5ng protein per band

Enabling

Mass Spec compatible
Edman Sequencing

Safe

Non toxic
No Methanol
Sink disposal

kem-en-tec-nordic.com

**KEM
EN
TEC
Nordic**



Does your gel stain tick all the boxes?

	InstantBlue™ (kem-en-tec)	SimplyBlue™ (Invitrogen)	GelCode™ Blue (Pierce)	Home-made (Anyone)
No need to wash before staining	✓	✗	✗	✗
No need to fix before staining	✓	✓	✗ ³	✓
No need to heat/microwave	✓	✗ ¹	✓ ⁴	✓ ¹
No need to de-stain	✓	✗	✗	✗
No risk of over-staining if left too long	✓	✗ ²	✗	✗
Non toxic and easily disposed of	✓	✓	✗	✗
True single step staining	✓	✗	✗	✗

¹ Recommended, otherwise requires 1 or more hours to stain ² Sensitivity decreases after 3 hours staining ³ If gels run with MOPS or MES buffer systems ⁴ Suggested to provide faster, 30 minute staining but sensitivity reduced

InstantBlue™ Ultrafast protein stain

For further information go to
www.kem-en-tec-nordic.com



Kem-En-Tec Nordic A/S
Kuldysen 10
DK-2630 Taastrup
Phone: +45 3927 1777

info@kem-en-tec.com
www.kem-en-tec-nordic.com